



Suzhou ColdStones Technology Co.,Ltd

Partial product summary

Name	Model	Diameter	Length	Purity	Appearance	Optional Solvent	Optional Concentration
silver nanowires	CST-NW-S40	40±10nm	35±5um	>99.5%	Gray suspension (Dispersed in solvent)	Water/Ethanol/IPA/Water plus binder(w+b)/Customized	5mg/ml、10mg/ml
	CST-NW-S70	70±10nm	35±5um				
	CST-NW-S90	90±10nm	40±5um				
	CST-NW-S120	120±20nm	40±5um				
	CST-NW-S250	250±50nm	40±5um				
	CST-NW-S400	400±100nm	40±5um				
Name	Model	Diameter	/	Purity	Appearance	Optional Solvent	Optional Concentration
silver nanoparticles	CST-NP-S2	2±0.5nm	/	>99.5%	Black suspension	Water/Ethanol/IPA/Water plus binder(w+b)/Customized	5mg/ml、10mg/ml
	CST-NP-S4	4±1nm	/		Blue suspension		
	CST-NP-S7	7±2nm	/		Blue suspension		
	CST-NP-S10	10±2nm	/		Light Yellow suspension		
	CST-NP-S20	20±3nm	/		Light Yellow suspension		
	CST-NP-S50	50±5nm	/		Black suspension		
	CST-NP-S80	80±5nm	/		Yellow suspension		
	CST-NP-S100	100±10nm	/		Grayish Yellow suspension		
	CST-NP-S150	150±10nm	/		Black suspension		
	CST-NP-S200	200±10nm	/		Grayish Yellow suspension		
Name	Model	Diameter	Thickness	Purity	Appearance	Optional Solvent	Optional Concentration
silver nanoflakes	CST-NF-S20	20±5nm	10-80nm	>99.5%	Light blue suspension	Water/Ethanol/IPA/Water plus binder(w+b)/Customized	5mg/ml、10mg/ml
	CST-NF-S40	40±10nm	10-80nm		Light blue suspension		
	CST-NF-S80	80±10nm	10-80nm		Aubergine suspension		
	CST-NF-S100	100±10nm	10-80nm		Pink suspension		
	CST-NF-S150	150±10nm	10-80nm		Reddish suspension		
	CST-NF-S400	400±50nm	10-80nm		Yellow suspension		
Name	Model	Average Diameter	Length	Purity	Appearance	Optional Solvent	Optional Concentration
copper nanowires	CST-NW-C20	20nm	30±5um	>99.5	brick-red suspension	Ethanol/IPA/N-hexane/Customized	5mg/ml、10mg/ml
	CST-NW-C60	60nm	30±5um		brick-red suspension		